

HER2 Discordant Results in Local vs. Central Testing in the Phase 3 Nelipepimut-S Trial and Implementation of Leica Bond Oracle HER2 Immunohistochemistry (IHC) System for Low and Intermediate Levels (1+, 2+) of HER2 Protein Expression as a Companion Diagnostic

M Melisko MD¹, E Mittendorf MD PhD², S Safina MD PhD³, M Schenker MD PhD⁴, AM Brunt MD⁵, M Litwiniuk MD⁶, J Mackey MD⁷, K Petrakova MD⁸, S Alieva MD⁹, L Chance¹⁰, GS Choy PharmD¹⁰, L Morgenthien PhD¹¹, S Kumar PhD¹², HS Rugo MD¹

¹University of California, San Francisco, San Francisco, CA, USA; ²University of Texas, M.D. Anderson Cancer Center, Houston TX, USA; ³Republican Clinical Oncology Center, Kazan, Republic of Tatarstan, Russian Federation; ⁴Sf Nectarie Oncology Center, Craiova, Romania; ⁵University Hospital of North Staffordshire, Stoke on Trent, UK; ⁶Greater Poland Cancer Center, Poznan, Poland; ⁷Cross Cancer Institute, Edmonton, Alberta Canada; ⁸Masaryk Memorial Cancer Institute, Cancer Care Clinic, Brno, Czech Republic; ⁹Donetsk Regional Antitumor Center, Donetsk, Ukraine; ¹⁰Galena Biopharma, Inc., Portland OR, USA; ¹¹ICON, NC, USA; ¹²Leica Biosystems, Danvers Massachusetts, USA.


Background

- ❖ Significant discordance exists between local and central laboratory testing for HER2 overexpression
- ❖ Discordance rates for tumors selected for HER2 0-2+ are not well described
- ❖ It is important to document sources and rates of discordance to avoid under and over treatment with HER2-targeted therapies

PRESENT Study

- ❖ NeuVax™ is a peptide vaccine that targets cells expressing any level of HER2
- ❖ The PH3-01, PRESENT (Prevention of Recurrence in Early Stage, Node Positive Breast Cancer with Low to Intermediate HER2 Expression with NeuVax™ Treatment) study is an international, multi-center randomized double-blind trial of NeuVax™ in 700 patients (NCT01479244)
- ❖ Patients are randomized to receive NeuVax™ + adjuvant vs. placebo + adjuvant over 36 months
- ❖ Study is open for enrollment in the US, Canada, UK, Germany, France, Israel, Russia, Ukraine, Poland, Romania, Hungary, Czech Republic, and Bulgaria
- ❖ **Objectives and Endpoints:**
 - Assess efficacy/safety of vaccine NeuVax™
 - Endpoints: DFS, OS, time to recurrence
- ❖ **Eligibility:**
 - Early-stage, node-positive breast cancer
 - Completed adjuvant chemotherapy within 12 wks
 - HLA-A2 or A3 positive
 - Low/intermediate expression of HER2 protein (1+ or 2+ by IHC; FISH non-amplified)
- ❖ **Methods: Screening for HER2 Expression:**
 - Local HER2 testing: samples for any degree of HER2 by local results to be assessed by central laboratory for confirmation of eligibility
 - Initial HER2 central testing: central laboratory confirmation using DAKO HercepTest™
 - Current HER2 central testing: upon approval of protocol amendment, patients are evaluated with Leica BOND Oracle III for HER2 eligibility

Figure 1. Comparison of Local vs. Central IHC (N = 1697)*



■ = No discordance or discordance did not affect eligibility
■ = Discordance resulted in change in eligibility
■ = Change in eligibility was to/from HER2 3+

¹**Eligibility definitions:** Local testing: Any degree of HER2 testing by IHC is eligible for screening
 Central testing: HER2 1+ or 2+ by DAKO HercepTest™ are eligible for enrollment

- ❖ Overall discordance rate = 54% (921/1697)
- ❖ Discordance resulted in changes in eligibility for 38% of patients
 - 311 patients (18%) were not enrolled due to discordance between local and central results

Table 1. Global IHC Discordance Rates (N = 1697)*

Local Testing	Central Testing (DAKO HercepTest™)				Discordance Rate
	0 (n = 482)	1+ (n = 731)	2+ (n = 307)	3+ (n = 177)	
0 (n = 584)	245 (42%)	247 (42%)	77 (13%)	15 (3%)	58%
1+ (n = 707)	182 (26%)	362 (51%)	132 (19%)	31 (4%)	49%
2+ (n = 300)	48 (16%)	114 (38%)	88 (29%)	50 (17%)	71%
3+ (n = 106)	7 (7%)	8 (8%)	10 (9%)	81 (76%)	24%
Discordance Rate	49%	50%	71%	54%	

Results

Table 2. IHC Discordance Rates by Location

Region	Country	Discordance Rate (n/N)	
North America	US	46% (86/185)	47% (102/219)
	Canada	47% (16/34)	
Eastern Europe	Bulgaria	70% (62/89)	55% (726/1309)
	Czech Republic	50% (55/109)	
	Hungary	58% (48/83)	
	Poland	60% (25/42)	
	Romania	59% (48/82)	
	Russia	53% (287/539)	
Western Europe	Ukraine	55% (201/365)	55% (79/143)
	France	38% (5/13)	
	Germany	60% (28/47)	
Middle East	UK	55% (46/83)	54% (14/26)
	Israel	54% (14/26)	

- ❖ Discordance (local vs. central) rates varied by geographic location (Table 2)
 - Ranged from 38% to 70% by country
 - Discordant rates were comparable and similar across regions
 - Rates of patients who were locally ineligible but proved to be eligible by central lab ranged from 29% (Israel) to 67% (UK, Canada)

LEICA Bond Oracle HER2 IHC Assay

- ❖ The Leica Bond Oracle HER2 IHC System is an FDA-approved fully-automated assay to accurately determine HER2 oncoprotein status in breast cancer tissue as an aid in the assessment of patients for whom Herceptin® treatment is being considered
- ❖ The Leica assay include four cell line controls, each one with a unique staining pattern representing the four HER2 classification (0, 1+, 2+, 3+)
- ❖ The Leica Bond Oracle HER2 IHC assay is also validated to reliably identify IHC HER2 1+ and 2+ cases in order to support the on-going PRESENT trial as a companion diagnostic for NeuVax™; whereas, the DAKO HercepTest is not validated to distinguish between HER2 IHC 0 and HER2 1+

Table 3. IHC DAKO vs. LEICA**

DAKO HercepTest™		Leica Bond Oracle			
		0 (n = 26)	1+ (n = 138)	2+ (n = 24)	3+ (n = 1)
DAKO HercepTest™	0 (n = 16)	6 (23%)	10 (7%)	0 (0%)	0 (0%)
	1+ (n = 117)	17 (65%)	91 (66%)	9 (38%)	0 (0%)
	2+ (n = 56)	3 (12%)	37 (27%)	15 (63%)	1 (100%)
	3+ (n = 0)	0 (0%)	0 (0%)	0 (0%)	0 (0%)

** Data cut-off: October 29, 2014

Conclusions

- ❖ Significant discordance exists between local and central laboratory test results for HER2 expression by IHC, even at lower levels of expression
- ❖ Discordance is similar across geographic regions
- ❖ Current available tests (DAKO HercepTest™) are defined by their ability to determine HER2 3+ status by IHC
- ❖ To improve accuracy and specificity for HER2 1+ and 2+ status testing and investigate a companion diagnostic for NeuVax™, the Leica Bond Oracle HER2 IHC System has been incorporated as central HER2 IHC screening for the PRESENT study
- ❖ Preliminary early LEICA results demonstrate high to moderate one-way concordance in HER2 1+/2+ compared to DAKO HercepTest™

The authors would like to thank the investigators and patients participating in the PRESENT study

The PRESENT study is sponsored by Galena Biopharma, Inc. This presentation is the intellectual property of the author/presenter. Contact at lchance@galenabiopharma.com for permission to reprint and/or distribute.

Disclosures:
 HS Rugo serves on the Scientific Advisory Board for Galena Biopharma, Inc. All other investigators receive research funding from Galena Biopharma, Inc. for the conduct of the PRESENT study.
 L Chance and G Choy are employed by Galena Biopharma, Inc.; L Morgenthien is employed by ICON.

